

P3 filter for CM-6 full face mask

COD ARTICOL

SM0704017999999



DESCRIZIONE

Possibility of filter sterilization up to 5 times, either soak in 60% alcohol, or by 80°C temperature for 20 minutes, or by ozone, or by UV. Nanofilters capture up to 50 nanometers, while the smallest virus has 200 nanometers. Much longer lifespan, thanks to the nanomembrane and one inhaling and exhaling valve around 4 times longer than usual half masks and full face masks filters.

STANDARDI

EN 143:2000/A1:2
006

P3R

P3 filter for CM-6 full face mask

CARACTERISTICI

Fabricat in EU
Greutatea articolului: 40g
Fără latex
Fără silicon

COD ARTICOL

SM0704017999999

STANDARDE

EN 143:2000/A1:2006 / P3R

Copyright © CERVA GROUP a.s .. Toate drepturile rezervate. Orice erori de imprimare, modificări de imagine și specificații ale produsului sunt rezervate.

Data modificării: 27 feb. 2023,
16:36:05

CERVA GROUP a.s.,
Průmyslová 483,
252 61 Jeneč,
Republica Cehă